




PRO

MAKE SURE YOUR AEDs WORK WHEN THEY ARE NEEDED

The Cardilink Connected Solutions PRO is a LIVE, 24/7, advanced active monitoring system. The device continuously checks the defibrillator to ensure it is fully functional.

It enables full integration, including deployment and motion detection, via a secure live IOT connection.



Product highlights

AED Readiness

Readiness indications to registered points of contact so that device owners have the confidence to deploy a device in an emergency or immediately address the needs of a device no longer ready.

AED Battery Status

Incorporates digital data readout providing exact battery percentage of the AED.

Motion Detection

Incorporates an accelerometer which detects the movement of the AED 24/7 ensuring that the AED is at its intended location.

Geo-Location

Features Geo-Location so that device owners may identify the whereabouts of deployed devices and ensure proper placement/location. This feature helps locate lost/misplaced devices or devices pending return from the SCA victim receiving hospital.

GPS-Tracking

Geo-tracking so that public devices may be tracked when moving away from the designated location (i.e - device theft and tracking of device movement).

4 Year Battery Life